



Ethernet based



Mounting pattern



4x screw M4x16 and 2x O-ring AS112 are part of the instrument delivery.

Seal surface
without cross scratch

Ra 0.8

Dimensions subject to change without notice.

Tolerances ± 0.2 unless
otherwise mentioned

Title: Dimensional drawing
FLEXI-FLOW Compact FF-a1nD


Bronkhorst®
WWW.BRONKHORST.COM

Projection:
Third Angle
Dimensions:
Metric (mm)

Drawn by: M.H.
Date: 21-09-2023
Rev. by:
Date:
Checked: V.H.
Date: 21-09-2023

Revision Description:
New drawing

Drawing no.: 7.15.285

Rev: A

Scale:
1 : 1
A3
Page:
1 - 1